

44070





Description

"Device functions: Strengthening muscles of upper parts, thighs abdomen and spine. It helps to maintain the correct posture. Preventive effect on the curvature of the spine. Method of use: Stand back to the device and place your hands on support bar, grab the handles. Pull up the legs To the body and then slowly leave."

Model typeSolitairSerial number44070Equipment nameLegs Lift

Join the action!



Legs Lift

373	Product line	Action Outdoor Fitness Solitair
73	Age range	10 - 99
	Unit measurements	62cm x 110cm x 205cm
7	Required surface area	375cm x 255cm
254	Impact area	10m²
	Free height of fall	0cm
	Material	Powder coated zinc layer steel. Pole is made of a steel tubes diameter 114 mm. Holders and other tube elements are made of steel tube diameter 40 mm, 3,2 mm thickness. Seats: weather resistant HDPE. Pedals: stainless steel.
	Heaviest part	0kg
	Largest part	0cm - 0cm
	Spare parts	Available
	Installation time without grondwork	2 hour(s) / 2 person(s)
	Special tools	Included
	Anchoring	Equipment will be installed 30 cm under ground level on a concrete foundation. Extra concrete is needed for more stability.
30	Remarks	Pylon (basic element) can be used for second device.
		l