



Description


"Device functions: It helps to build muscles of the lower limbs and waist. Affects the flexibility of the joints and improve blood circulation.
Method of use: Sit on the seat, rest your feet on the pedals. Straighten the legs pushing body away from the device and bend them at knees again."

Model type
Serial number
Equipment name

Solitaire
44100
Legs press

Join the action!

Legs press



| | |
|--|--|
| Product line | Action Outdoor Fitness Solitair |
| Age range | 10 - 99 |
| Unit measurements | 62cm x 115cm x 205cm |
| Required surface area | 360cm x 310cm |
| Impact area | 12m ² |
| Free height of fall | 0cm |
| Material | Powder coated zinc layer steel. Pole is made of a steel tubes diameter 114 mm. Holders and other tube elements are made of steel tube diameter 40 mm, 3,2 mm thickness. Seats: weather resistant HDPE. Pedals: stainless steel. |
| Heaviest part | 0kg |
| Largest part | 0cm - 0cm |
| Spare parts | None |
| Installation time without grondwork | 2 hour(s) / 2 person(s) |
| Special tools | None |
| Anchoring | Equipment will be installed 30 cm under ground level on a concrete foundation. Extra concrete is needed for more stability. |
| Remarks | Pylon (basic element) can be used for second device. |

